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# Preparation Instructions for the SDS1.2/ SDS2.2 Implant System

## Application:

Reusable instruments of the *SDS Instrument Set Implantology*, *SDS surgical tray* and single-use prosthetic components of the SDS implant system.

All instruments, the *SDS surgical tray* and the prosthetic components mentioned in these instructions are supplied non-sterile. They must be thoroughly cleaned, disinfected and sterilized before each use. None of the instruments or prosthetic components must be sterilized in the SDS original packaging.

## 1. Storage and handling

Before initial processing, all system components must be stored at room temperature in their original packaging under clean, dry and dust-free conditions. Instruments must be protected against external influences like impact, shock and falling when transported in the facility.

## 2. Cleaning/ disinfection/ sterilization

All instruments, the *SDS surgical tray* and the prosthetic components mentioned in these instructions are supplied non-sterile. The instruments and the *SDS surgical tray* are reusable, the prosthetic components may only be used once on a single patient. All system components must be thoroughly cleaned, disinfected and sterilized prior to first and subsequent use.

The *SDS surgical tray* must be cleaned, disinfected and sterilized before the first application and after each application. Cleaning/ disinfecting the equipped *SDS surgical tray* is not permitted.

Instruments and prosthetic components made of zirconium oxide must always be prepared separately from metal instruments, otherwise metallic abrasion may occur on the zirconium oxide surface. Therefore, only use brushes with metal-free bristles for pre-cleaning.

## Initial disinfection

1. Ensure that all instruments contaminated during surgery are collected separately. Do not return them to the *SDS surgical tray* to avoid contamination of the filled tray.
2. Immediately after use, all instruments must be immersed in a bath with combined cleaning and disinfecting agent. This serves for user safety and prevents the contaminants (blood, secretion, tissue remains) from drying.
3. The instruments are to be disinfected with disinfectants and cleansing agents with proven efficacy for disinfection (e.g. VAH/ DGHM-listed, CE-marking for EU or EPA-registered and FDA-compliant for USA), suitable for cleaning and disinfecting dental instruments and compatible with the materials of the products to be cleaned and disinfected. The manufacturer's instructions for the use of those disinfectants and cleansing agents, as for example concentrations and application times as well as instructions for post-rinsing have to be adhered to.

## Pretreatment after initial disinfection

1. Disassemble the *SDS insertion tools* by unscrewing the fastening screw counterclockwise. Disassemble the *SDS torque ratchet* according to the instructions for use.
2. Remove the instrument carrier from the *SDS surgical tray* and remove any instruments remaining in the instrument carrier.
3. Remove all contamination from the instruments/ instrument carrier/ *SDS surgical tray* under running water using brushes with metal-free bristles.
4. Full rinse of the instruments for at least 1 minute with deionized water with a low bacterial- (max. 10 bacteria/ ml) and low endotoxin count (max. 0.25 endotoxin units/ ml).
5. Additionally, rinse the lumen of the *SDS insertion tools* with deionized water with a low bacterial- (max. 10 bacteria/ml) and a low endotoxin count (max. 0.25 endotoxin units/ ml) at least five times using a syringe (at least 10 ml each).

## As part of your responsibility for the sterility of the products in application, please observe the following:

- in general, use only adequately validated methods specific to the equipment and product for cleaning/ disinfection and sterilization
- use only thermal disinfection that requires no disinfecting agent. There is a risk of disinfectant residues on the products when chemical disinfection is used
- regularly check and service the equipment used (thermal disinfectant, sterilizer)
- observe the validated parameters in each cycle

SDS recommend according to the recommendation of the Robert-Koch Institut the mechanical cleaning and disinfection by a standard automatic cleaning program in a washer-/disinfectant with proven efficacy (CE marking, compliant with EN ISO 15883-2 in Europe or AAMI TIR30:2011 (R) 2016 in the U.S.) using a proven program for thermal disinfection (AO value > 3.000 or - for older devices - at least 5 minutes at 90 °C/ 194 °F) which is suitable for the system components of the SDS implant system and includes sufficient rinsing cycles with deionized water with a low bacterial- (max. 10 bacteria/ ml) and low endotoxin count (max. 0.25 endotoxin units/ ml).

## Mechanical cleaning and disinfection

1. Place instruments/ emptied instrument carrier/ *SDS surgical tray* and prosthetic components in the washer/ disinfectant in such a way that the products are directly hit by the spray jet and don't come in contact to each other.
2. Put chemical detergent into the washer/ disinfectant, following the instructions of the manufacturer of the washer/ disinfectant.
3. Start the Vario TD program including thermal disinfection (see general note below for details of the program used for validation).
4. After program end remove instruments/ emptied instrument carrier/ *SDS surgical tray* and prosthetic components from the washer/ disinfectant and dry (preferably with filtered compressed air as recommended by the RKI).
5. Visual examination to ensure that the products are clean, undamaged and not corroded. In case of residual contamination after mechanical reprocessing, repeat the cleaning and disinfecting process until no visible contamination is left. In case of damage or corrosion, remove products from use.

## Sterilization in the autoclave: All non-sterile packed products must not be sterilized in the original packaging!

### Packaging:

1. Mount the *SDS insertion tools* by screwing in the fixing screw clockwise. Maintain and mount the *SDS torque ratchet* according to the instructions for use.
2. SDS recommend to place all components of the *SDS Instrument Set Implantology* after cleaning, disinfection and inspection in the *SDS Surgical Tray* and pack the tray with the

- instruments (or the individually packed instruments).
3. Pack the prosthetic components individually.
4. We recommend using disposable sterilization packaging with CE marking in Europe or FDA clearance in the U.S.
5. Ensure that the sterilization packaging is suitable for steam sterilization (constant temperature of at least 141°C/ 286°F, sufficient vapor permeability) and that the system components are adequately protected against mechanical damage.
6. When sealing the system components in the foil, make sure that the packaging is large enough to ensure that there is no pressure on the seal.

## Sterilization:

For this procedure, the following instructions have to be observed: Steam sterilization using validated fractionated vacuum processes according to ISO 17665 in a device that complies with the provisions of DIN EN 13060 or EN 285 in Europe or have FDA clearance in the U.S.

- Fractionated pre-vacuum (type B)
- Sterilization temperature: EU: 134 °C (273 °F) / U.S.: 132 °C (270 °F)
- Hold time: at least 4 minutes (full cycle)
- Drying time: at least 20 minutes

In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by the standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several products.

Follow the instructions of the device manufacturer. System components have to be checked for superficial damages (formation of cracks) and corrosion after sterilization. Do not reuse damaged system components.

The operator of medical products is responsible for making sure that cleaning, disinfection and sterilization processes are carried out by qualified personnel, using the appropriate materials and suited equipment.

**For the U.S.:** the validated procedures require the use of FDA-cleared sterilizers, sterilization trays, sterilization wraps, biological indicators, chemical indicators and other sterilization accessories labeled for the sterilization cycle recommended. The health care facility should monitor the sterilizer according to an FDA recognized sterility assurance standard such as ANSI/ AAMI ST79:2017/ AAMI TIR12:2010. Instruments must be sterilized.

**General note:** The proof of general suitability for effective mechanical cleaning and disinfection and sterilization was provided by an independent accredited testing laboratory taking into account the above-described procedure.

In detail, the validated cleaning and disinfecting process uses the thermal disinfectant (cleaning and disinfecting device) G 7835 CD (Miele & Cie. GmbH & Co., Gütersloh) and neodisher® MediClean forte (0.5 %) as the cleaning agent (Dr. Weigert GmbH & Co. KG, Hamburg) with parameters described below:

- 60 seconds pre-rinse (cold water of potable quality)
- 180 seconds pre-rinse (cold water of potable quality)
- 600 seconds clean (55 °C/ 131 °F water of potable quality and cleaner neodisher® MediClean forte (0.5%) (Dr. Weigert GmbH & Co. KG, Hamburg)
- 180 seconds rinse (cold deionized water)
- 300 seconds thermal disinfection (AO value > 3.000 or minimum 90 °C/ 194 °F)
- 30 minutes drying (min. 100 °C/ 212 °F)

The validated sterilization process refers to the universal program of the Tuttnauer company in the autoclave Tuttnauer EHS 3870 (fractionated vacuum process with 3 pre-vacuum phases, half-cycle hold time (132 °C/ 270 °F), drying time at least 20 minutes). The procedure described above was taken into account.

## 3. Storage

Store sterilized instruments according to the instructions of sterile packaging manufacturer.

## 4. Service life and reprocessing cycles

The service life of the instruments of the *SDS Instrument Set Implantology* depends on the bone quality and individual work habits of the user. 50 application cycles should not be exceeded. Prosthetic components may only be used once for a single patient.

## 5. Safety and liability:

Worn or damaged instruments or system components are to be discarded and immediately replaced by new ones. The above instructions are to be strictly observed.

**For technical support and further information please contact:**  
 SDS Swiss Dental Solutions AG  
 Konstanzerstrasse 11  
 CH-8280 Kreuzlingen