**Instrumentenset Implantologie**  
314223 - Rev 2021-09

Legende: Symbole der Etiketten · Legend: Symbols on the labels

**REF** Bestellnummer  
Catalog number

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## Instructions for the use of SDS Instrument Set Implantology

Article number/application:

Reusable dental drills for the preparation of drill holes in jawbone to insert dental implants.

### ⚠ Warning:

All instruments are delivered non-sterile and must be thoroughly cleansed, disinfected and sterilized before first and every use.

#### 1. Indication

The *SDS Instrument Set Implantology* is intended for the preparation and enlargement of the implant site. When using a drilling template, make sure to prevent the instrument from jamming to avoid possible fracture of the instrument.

#### 2. Contraindication

The patient must be informed of risks, side-effects and complications, as well as of necessary precautions in connection with *SDS Instrument Set Implantology*.

Known risk factors:

- Danger to anatomic structures in the region of the intended intervention
- Overheating of bone
- Drill fracture
- Fracture of instruments caused by excessive force/misuse by user

Side effects/possible adverse reactions:

- Pain, swelling, infection of soft- and hard tissue
- Dysaesthesia/paraesthesia
- No osseointegration
- Bone defects necessitating bone grafting
- Perforation of the sinus, mandibular base, floor of the mouth or lower alveolar ducts
- Damage to neighboring teeth/tooth roots
- Excessive bone loss, which might necessitate surgical intervention
- Aesthetic problems

#### 3. Treatment sequence/appropriate use

The product is intended for use by qualified dentists. General instructions on various surgical techniques can be found in the relevant specialist literature. Prior to the surgical intervention, patients have to be informed on the generally valid precautionary measures and rules of conduct. In addition to the orthopantomogram, we recommend a computer tomographic examination to determine the exact position and depth of the drill hole. Examine the adjacent structures to exclude any danger to adjacent structures. Administer a local anesthesia to the site where the instrument is to be used. Cut and prepare the mucosa and the periosteum according to the surgical procedure for implants. The exact position of the alveolar cavity is determined by means of the round bur (*SDSrb230*). The pilot drill (*SDSpd250*) is applied at the site previously determined by the round bur. The pilot drill serves to drill the initial hole for the intended preparation. When using drilling templates, make sure to prevent the instrument from jamming to avoid possible fracture of the instrument. Form- and compact drills (*SDSsd300Ord*, *SDSsd370rd*, *SDSsd470rd*, *SDSsd300dt*, *SDSsd380dt*, *SDSsd460dt*, *SDSsd350C*, *SDSsd430C*, *SDSsd510C*, *SDSsd670C*) are used for enlarging the cavity up to its final width. Drill in an intermittent mode, providing constant external cooling with a precooled (5°C/41°F) physiological saline solution. The external irrigation prevents thermal damage to the bone tissue and facilitates the evacuation of bone chips. Make sure that the drill does not jam or get stuck during use (increased risk of fracture). Apply low pressure during the preparation up to the required depth. Work at a speed of 300–600 rpm. Do not exceed the recommended speed to prevent fracture of the instrument. In the cortical region, the cavity can be enlarged with the countersink (*SDScs500*, *SDScs600*, *SDS2.0\_CS-500*, *SDS2.0\_CS-600*, *SDS2.0\_CS-600+*). SDS drilling sequences must be followed (see separate SDS user manuals).

#### 4. Accessories

Instruments with contra-angle shank can be connected to the shank extension (*SDSse001*) which extends the instrument by 15 mm. To this end, the locking hook of the shank extension must click into place with an audible sound. The drilling depth and direction can be checked with the help of the direction and depth indicator (*SDSdg240*). The instrument is placed in the alveolar cavity and the depth is determined by means of the ring-shaped depth marks.

#### 5. Storage and handling

Before initial processing, all system components must be stored at room temperature in their original packaging under clean, dry and dust-free conditions. Store sterilized instruments according to the instructions of sterile packaging manufacturer, preferably in the *SDS surgical tray*. Instruments must be protected against external influences like impact, shock and falling when transported in the facility.

#### 6. Cleaning/disinfection/sterilization

All instruments mentioned in these instructions are reusable and supplied non-sterile. They therefore have to be thoroughly cleansed, disinfected and sterilized prior to first- and every consecutive use. All instruments must not be sterilized in the SDS original packaging!

Instruments made from zirconia always must be processed separately from metal instruments, otherwise metallic wear on the zirconia instrument surface could originate. Only use brushes with metal-free bristles for pre-cleaning ceramic instruments as these do not leave discolorations caused by abrasion.

#### Initial disinfection

1. Immediately after use, all instruments must be immersed in a bath with combined cleaning and disinfecting agent. This serves for user safety and prevents the contaminants (blood, secretion, tissue remains) from drying.
2. The instruments are to be disinfected with disinfectants and cleansing agents with proven efficacy for disinfection (e.g. VAH/DGHM-listed, CE-marking for EU or EPA-registered and FDA-compliant for USA), suitable for cleaning and disinfecting dental instruments and compatible with the materials of the products to be cleaned and disinfected. The manufacturer's instructions for the use of those disinfectants and cleansing agents, as for example concentrations and application times as well as instructions for post-rinsing have to be adhered to.

#### Pretreatment after initial disinfection

1. Remove all contamination from the instruments under running water using brushes with metal-free bristles.
2. Full rinse of the instruments for at least 1 minute with deionized water with a low bacterial- (max. 10 bacteria/ml) and low endotoxin count (max. 0.25 endotoxin units/ml).

As part of your responsibility for the sterility of the products in application, please observe the following:

- in general, use only adequately validated methods specific to the equipment and product for cleaning/disinfection and sterilization
- use only thermal disinfection that requires no disinfecting agent. There is a risk of disinfectant residues on the products when chemical disinfection is used
- regularly check and service the equipment used (thermal disinfectant, sterilizer)
- observe the validated parameters in each cycle

SDS recommend according to the recommendation of the Robert-Koch Institut the mechanical cleaning and disinfection by a standard automatic cleaning program in a washer/-disinfectant with proven efficacy (CE marking, compliant with ISO 15883-1 and ISO 15883-2 in Europe or ANSI/AAMI ST 15883-1/AAMI TIR30 in the U.S.) using a proven program for thermal disinfection (A<sub>0</sub> value > 3.000 or – for older devices – at least 5 minutes at 90°C/194°F) which is suitable for the *SDS Instrument Set Implantology* and includes sufficient rinsing cycles with deionized water with a low bacterial- (max. 10 bacteria/ ml) and low endotoxin count (max. 0.25 endotoxin units/ml).

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#### Mechanical cleaning and disinfection

1. Place instruments in the washer/disinfectant in such a way that the products are directly hit by the spray jet and don't come in contact to each other.
2. Put chemical detergent into the washer/disinfectant, following the instructions of the manufacturer of the washer/disinfectant.
3. Start the automatic cleaning program including thermal disinfection (see general note below for details of the program used for validation).
4. After program end remove instruments from the washer/disinfectant and dry if necessary (preferably with filtered compressed air as recommended by the RKI).
5. Visual examination to ensure that the products are clean, undamaged and not corroded. In case of residual contamination after mechanical reprocessing, repeat the cleaning and disinfecting process until no visible contamination is left. In case of damage or corrosion, remove products from use.

Sterilization in the autoclave: ⚠ All non-sterile packed products must not be sterilized in the original packaging!

#### Packaging

1. SDS recommend to place all components of the *SDS Instrument Set Implantology* after cleaning, disinfection and inspection in the *SDS Surgical Tray* and pack the tray with the instruments (or the individually packed instruments).
2. Use prefabricated sterile barrier system according to the requirements of ISO 11607-1/ANSI/AAMI ISO 11607-1 (disposable sterilization packaging with CE marking in Europe or FDA clearance in the U.S)
3. Ensure that the sterilization packaging is suitable for steam sterilization (constant temperature of at least 141°C/286°F, sufficient vapor permeability) and that the instruments are adequately protected against mechanical damage.
4. When sealing the instruments in the foil, make sure that the packaging is large enough to ensure that there is no pressure on the seal.

#### Sterilization

For this procedure, the following instructions have to be observed: Steam sterilization using validated fractionated vacuum processes according to ISO 17665-1/ANSI/AAMI ISO 17665-1/ANSI/AAMI/ST79 in a device that complies with the provisions of DIN EN 13060 or EN 285 in Europe or have FDA clearance in the U.S. Maximum sterilization temperature may not exceed 138°C (280°F).

- Fractionated pre-vacuum (type B)
- Sterilization temperature: EU: 134°C (273°F)/U.S.: 132°C (270°F)
- Hold time: at least 4 minutes (full cycle)
- Drying time: at least 20 minutes

In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by the standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several products. Follow the instructions of the device manufacturer.

Ceramic instruments have to be checked for superficial damages (formation of cracks) after sterilization. Do not reuse damaged instruments due to the increased risk of fracture.

The operator of medical products is responsible for making sure that cleaning, disinfection and sterilization processes are carried out by qualified personnel, using the appropriate materials and suited equipment.

For the U.S: the validated procedures require the use of FDA-cleared sterilizers, sterilization trays, sterilization wraps, biological indicators, chemical indicators and other sterilization accessories labeled for the sterilization cycle recommended. The health care facility should monitor the sterilizer according to an FDA recognized sterility assurance standard such as ANSI/AAMI ST79:2017. Instruments must be sterilized.

General note: The proof of general suitability for effective mechanical cleaning and disinfection and sterilization was provided by an independent accredited testing laboratory taking into account the above-described procedure.

In detail, the validated cleaning and disinfecting process uses the thermal disinfectant (cleaning and disinfecting device) G 7835 CD (Miele & Cie. GmbH & Co., Gütersloh) and neodisher® MediClean forte (0.5 %) as the cleaning agent (Dr. Weigert GmbH & Co. KG, Hamburg) with parameters described below:

- 60 seconds pre-rinse (cold water of potable quality)
- 180 seconds pre-rinse (cold water of potable quality)
- 600 seconds clean (55°C/131°F water of potable quality and cleaner neodisher® MediClean forte (0.5%) (Dr. Weigert GmbH & Co. KG, Hamburg))
- 180 seconds rinse (cold deionized water)
- 300 seconds thermal disinfection (A<sub>0</sub> value > 3.000 or minimum 90°C/194°F)
- 30 minutes drying (min. 100°C/212°F)

The validated sterilization process refers to the universal program of the Tuttnauer company in the autoclave Tuttnauer EHS 3870 (fractionated vacuum process with 3 pre-vacuum phases, half-cycle hold time (132°C/270°F), drying time at least 20 minutes). The procedure described above was taken into account.

#### 7. Service life and reprocessing cycles

Wear of drills is dependent on bone quality and individual work habits of the user. 50 application cycles should not be exceeded.

#### 8. Disposal

Adhere to the general requirements for the disposal of medical devices when disposing of SDS products or packaging material.

#### 9. Safety and liability:

Worn and damaged instruments are to be discarded and immediately replaced by new ones. The above instructions regarding the handling of the instruments are to be strictly observed. Only use the instruments according to their intended purpose. The non-observance of these safety instructions can lead to injuries.

It is the responsibility of the user to check the products prior to use to ensure that they are suited for the intended purpose. In case of contributory negligence by the user, SDS partially or totally declines liability for all resulting damages, particularly due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

Warranty: The SDS product may only be used according to the manufacturer's instructions. The operator is responsible for ensuring that the product is used for its intended purpose and must also assess whether the product is suited to the patient's particular situation. The products may only be used together with SDS products. The SDS warranty is invalidated by the use of third-party products that are not approved by SDS. Liability will not be accepted for products that have been modified, misused or fitted incorrectly.