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 MOJE KERAMIK-IMPLANTATE GmbH & Co. KG
Eschenweg 10-12
D-07616 Petersberg

 CE 0483

Instructions for Use for Two-Piece SDS2.1 Dental implant system

Description: Two-piece SDS2.1 implants are inserted into the jawbone for transgingival healing. The separately packaged SDS2.1 standard implant post is used for fixing the prosthetic restoration. For this purpose, the SDS2.1 standard implant post is cemented with glass ionomer cement (Ketac Cem™) after the implant interior geometry has been cleaned and dried.

The SDS2.1 standard implant post can be cemented either directly after implant placement or later after implant healing. SDS2.1 standard implant posts must be disinfected and sterilized before cementation.

Intended Purpose: The dental implants from the manufacturer are used as a replacement in case of destruction or loss of the tooth root. The implants are permanent implants designed to replace normal body structures. The cause of the destruction of the root of the tooth is irrelevant. Essentially, the implant is an artificial substitute for a natural root.

Packaging/Material: Sterile individual package: SDS2.1 dental implant is made of high strength zirconia TZP-A (tetragonal zirconia polycrystal) zirconium dioxide ceramics (bio ceramics under ISO 13356)

Unsterile individual package: SDS2.1 standard implant post is made of high strength zirconia TZP-A (tetragonal zirconia polycrystal) zirconium dioxide ceramics (bio ceramics under ISO 13356)

Unsterile individual package: silicone tubes made of poly-dimethylsiloxane

**⚠ The products are single-use disposable devices and not intended for reuse!
Unsterile delivered products must be reprocessed by an appropriate procedure before use (see section "Validated treatment process").**

**Indications SDS2.1:
Please check detailed clinical instructions in the SDS2.1 product information brochure.**

SDS2.1 dental implants are surgically inserted into the upper and lower jaw to replace missing teeth. They serve as attachment of prosthodontic appliances and are particularly suitable for patients with metal allergies and the chronic diseases resulting from them.

1. Indicated for situations where implants are connected by interlocking or bridge restoration
2. **SDS2.1_38XX implants** are approved only as single tooth implant for upper lateral and lower incisors, as well as for premolars and splinted implants
3. **SDS2.1_46XX implants** are approved as single-tooth implants for incisors, cuspids, premolars and molars and for bridge fittings
4. SDS2.1 dental implants must be inserted at the tissue-level

Contraindications SDS2.1 dental implant system:

Pre-existing medical conditions or poor general condition can limit the possibility to insert dental implants surgically. Bruxism and insufficient bone quality/-quantity require special measures to ensure treatment success. SDS dental implants are not suitable for indications where there is a risk of excessive flexural moment (e.g. bridges with more than one pontic crowns/bridges with cantilevers).

1. Implants with small diameter and angled standard implant post are not recommended for the posterior region

2. SDS2.1 dental implants are not approved for bone-level position
3. SDS2.1_38XX implants are not approved for upper central incisors, canines, molars
4. SDS2.1_38XX implants are not approved for bridge fittings
5. Bone not fully healed (residual osteitis/NICO)
6. Severe general health conditions in the patient
7. Bruxism
8. Untreated periodontitis, poor oral hygiene, untreated abscess or bone infection
9. Crown-length longer than osseointegrated thread segment
10. Cantilever bridges (mesial or distal)
11. Pontic width between two bridge abutments greater than the width of one premolar
12. Connection of a natural tooth with an implant

⚠ Warning: Anatomic and general medical conditions can have a negative effect on the performance of dental implants. If the implant is exposed to pressure beyond its capacity, excessive bone loss or fracture of the implant can occur.

The dentist/oral surgeon is obliged to inform patients about the following risk factors:

- Poor bone quality
- Poor oral hygiene
- Diseases such as blood diseases or untreated hormonal disorders
- Alcohol or drug abuse
- Two-piece SDS2.1 implants must be protected from load during healing

⚠ For laser applications, the dentist/oral surgeon must follow the instructions of the device manufacturer

Possible adverse reactions: The following adverse reactions are possible after insertion of a dental implant:

1. Pain, swelling, infection of soft and hard tissue
2. Dysaesthesia/paraesthesia
3. No osseointegration
4. Loss of osseointegration
5. Bone defects necessitating bone reconstruction
6. Perforation of the maxillary sinus, base of mandible, base, floor of the mouth or lower alveolar ducts
7. Damage to adjoining teeth/tooth roots
8. Excessive bone loss that may require surgical intervention
9. Esthetic problems
10. Implant fracture

⚠ Precautionary measures: Intensive diagnostics of the oral cavity must be performed ahead of every implantation. It is necessary to take and assess the bone structures (OPG/CBCT/CT) must be performed and analyzed. When selecting the appropriate implant size/type, expected physiological chewing forces, as well as any parafunctional habits, must be taken into consideration.

Special training by SDS is recommended to guarantee safe and effective insertion of ceramic SDS implants. Surgical standard procedures must be complied with, maximum insertion torque of 35 Ncm may not be exceeded. Incorrect surgical techniques can lead to functional failure of the implant and bone loss of the supporting bone structure or other side effects.

SDS2.1 implants may only be restored with the enclosed SDS2.1 standard implant posts. SDS implants may only be restored further with SDS products (or products which are

authorized by SDS).

The silicone tube is intended to fix the standard implant post during cementation. The tube is pulled over the standard implant post before cementing. By gripping the tube with dental forceps, the standard implant post can be easily inserted into the implant and cemented.

The tip of the tweezers is guided through the silicone tube and pressed onto the abutment during cementation in order to ensure that the standard implant post is fully inserted in the implant.

The silicone tube may only be removed after the cement is hardened completely (Note the manufacturer's specifications of the cement).

⚠ Please also refer to information on cementing in the manufacturer's instructions for use

Treatment process for products delivered non-sterile:

⚠ The products must not be sterilized in the original packaging!

SDS2.1 standard implant posts are supplied non-sterile and must be cleaned, disinfected and sterilized before use. Permissible steam sterilization methods are fractionated vacuum methods (with sufficient product drying). Other sterilization methods (including steam sterilization using the gravity method) are not permitted.

Validated treatment process:

Cleaning:

Cleaning and disinfecting device:

"G 7835 CD" washer-disinfector (Miele & Cie. GmbH & Co., Gütersloh)

neodisher® MediClean forte (0.5 %) cleaning agent (Dr. Weigert GmbH & Co. KG, Hamburg)

- Pre-rinse for 60 seconds (cold water of drinking water quality)
- Pre-rinse for 180 seconds (cold water of drinking water quality)
- Clean for 600 seconds (55 °C/131 °F, water of drinking water quality and cleaning agents mentioned above Rinse for 180 seconds (cold, deionized water)
- Thermal disinfection for 300 seconds (A0 > 3000 or at least 90 °C/194 °F)
- Dry for 30 minutes (at least 100 °C/212 °F)

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



Universal program: fractionated vacuum process with 3 pre-vacuum phases, half-cycle dwell time (132 °C/270 °F), minimum drying time 20 minutes

Storage and handling: SDS dental implants are supplied in sterile packaging and must be stored in the original packaging under clean and dry conditions. If moisture is present in the inner packaging or if the package is damaged or has been unintentionally opened, the implants must not be used. SDS dental implants must not be sterilized again. Before use, the size/type and expiration date of the implants must be checked against the label on the package. Appropriate hygienic measures must be observed when handling the implants.

Contact with objects which could damage the implants must be avoided. SDS will not accept return of goods in damaged packages.

Disposal: For the disposal of SDS implants, the packaging material and the accessories, the general requirements for disposal of medical devices must be considered.

GRAPHIC SYMBOLS:

	CATALOG NUMBER
	BATCH CODE
	UNIQUE DEVICE IDENTIFIER
	MEDICAL DEVICE
	CONSULT INSTRUCTIONS FOR USE
	STERILIZED USING STEAM
	DOUBLE STERILE BARRIER SYSTEM
	DO NOT RESTERILIZE
	DO NOT RE-USE
	USE-BY DATE
	MANUFACTURER
	DATE OF MANUFACTURE
	DISTRIBUTOR
	CE MARKING OF CONFORMITY
	CE MARKING OF CONFORMITY NB NUMBER 0483
	NON-STERILE
	CAUTION
	TEMPERATURE LIMIT 10°C TO 35°C
	KEEP DRY
	KEEP AWAY FROM SUNLIGHT
	DO NOT USE IF PACKAGE IS DAMAGED AND CONSULT INSTRUCTIONS FOR USE
	AUTHORIZED REPRESENTATIVE IN THE SWITZERLAND

The detailed directions for use in the product brochure must be complied with without fail and are a basic requirement for the treatment.

For technical support and further information please contact:

 CH REP

SDS Swiss Dental Solutions AG
Konstanzerstrasse 11
CH-8280 Kreuzlingen

 CE 0483

MOJE KERAMIK-IMPLANTATE GmbH & Co. KG
Eschenweg 10-12
D-07616 Petersberg