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English

Product overview

Torque wrench set, consisting of :
Torque wrench
Adjustment key
Lubricant
Instructions for use

Description of the torque wrench

The torque wrench with adjustable force is a dental device used to tighten or loosen screws, prosthetic components and implants. It is a precision instrument that can be disassembled and that is supplied unsterile. To ensure that it functions perfectly every time, the torque wrench must be disassembled, disinfected, cleaned, greased and sterilised after each use, according to the instructions for use. Persons who use this device are obliged to read the instructions for use prior to handling. The handling and the use of the product are carried out without direct control from our side and remain under responsibility of the user. The user is liable for any possible damage that could occur. Before each use, in order to guarantee high torque precision, the device must be checked for correct function. This instrument is not a measuring device.

Use

By turning the torque adjustment screw, the torque wrench can be set to the desired torque value. To set the torque value correctly, the torque adjustment screw must be turned clockwise to reach the required torque value and set to the exact line marking. Ensure that the line on the handle is in straight alignment with the line on the torque adjustment screw. In order to change from a higher to a lower torque value, one must screw two turns under the desired torque value, then screw clockwise again to the exact line marking.

Ratchet mode could be set by turning the torque adjustment screw until the R marking for the torque wrench including Ratchet option.

The word 'IN' on the cover (3), shows the position of the wrench that is used for tightening ; by turning the device over, the word 'OUT' is used for loosening screws.

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± 3,5 Ncm with total confidence of 95%

Lifetime

3 years after supply date of manufacturer

Recommendations

This instrument must not be used for any applications other than those listed in the section "Description of the torque wrench" or with equipment that could damage the intended use of the device.

The persons in charge for the use and maintenance of this dental instrument should monitor any deterioration of the tightening, ratchet and torque mechanism of the device and, in the event of a defect, return the wrench to the supplier.

During assembly, it is essential not to mix the various components belonging to different torque wrenches because the components are not interchangeable.

If a component is lost, please return the whole instrument immediately to your retailer for repair.

Components cannot be sold separately.

Do not store the wrench with the spring compressed but with the torque set to its minimum.

This device must not be sterilised in the packaging provided by the manufacturer

Validity

This document cancels and replaces all previous versions.
The specifications contained herein can change without warning.
The official version of this document is the version in French.

Individual parts of the torque wrench


Fig. 1 The torque wrench is composed as follows :

- | | |
|--------------------------|----------------------------|
| 1. ratchet wheel | 6. spring / stop |
| 2. ratchet | 7. torque adjustment screw |
| 3. cover | 8. hexagon shank bit |
| 4. ratchet head / handle | 9. bolster (PPS) |
| 5. screw | 10. brake (PTFE) |

Cleaning the torque wrench

When used in situations that do lead to operative residues (blood, secretions, tissue remnants), the torque wrench must be disassembled completely (fig. 2 to 4) and placed in a suitable bath of disinfection in accordance with the recommendations of the manufacturer.

This operation facilitates cleaning because dry residues cause corrosion. After cleaning, thoroughly rinse the parts with water and use a nylon brush to rub internal and external surfaces of the various parts of the torque wrench. During the cleaning process, avoid all contact between each part of the torque wrench.

Disassembling the torque wrench

Fig. 2
Slide spring out of the handle

Fig. 2
Unscrew completely the torque adjustment screw (7) and extract the spring/stop (6).

If necessary, use the hexagon shank bit (8) to extract the spring. Be careful not to damage the hexagon shank bit. Do not separate the stop from the spring (6).

Fig. 3
Unscrew the screw (5) with the hexagon shank bit

Fig. 3
With the hexagon shank bit (8) of the torque adjustment screw (7), unscrew the screw (5), while maintaining a light pressure, to remove the cover (3).

Fig. 4
Take out the ratchet and the ratchet wheel from the ratchet head.

Fig. 4
The ratchet head is composed of two removable elements: the ratchet (2) and the ratchet wheel (1). The screw (5) is fixed inside the ratchet head.

Assembling the torque wrench

Fig. 5
Insert the ratchet wheel and the ratchet into the ratchet head


Fig. 5
To assemble the torque wrench, insert the next two elements in the following sequence: first the ratchet wheel (1) and then the ratchet (2).

Fig. 6
Lubricating point (▲)


Fig. 6
Moderately lubricate the functional contact area between the teeth of the ratchet wheel (1) and the pivot point of the ratchet (2) (▲) as shown in fig. 6. Remove any trace of the lubricant from the outer surface of the torque wrench because excessive lubricant leaves traces on the instrument's surface during sterilisation. Use only the instrument lubricant delivered with the torque wrench.

Fig. 7
Screw the cover on the ratchet head

Fig. 7
Once the parts (1) and (2) are in place, put the cover (3) on the ratchet head (4) and adjust. Turn the torque wrench while holding the cover. Fix the screw (5) tightly with the hexagon shank bit (7) of the torque adjustment screw.

Fig. 8
Insert spring / stop inside the handle and screw the torque adjustment screw


Fig. 8
The assembly of the torque wrench is complete as soon as the spring/stop is in place and the torque adjustment screwed back in the handle. Check the function by activating the ratchet wheel. If there is a regular ratchet noise, the instrument operates correctly. Before use, always check the assembly and the correct function of the torque wrench.

Fig. 9
Use the adjustment key (provided) to obtain the higher torque values more easily

Sterilisation

Permissible steam sterilization procedures are fractionated vacuum procedures (with sufficient product drying). Other sterilization procedures (including steam gravity sterilization) are not allowed. Care should be taken that:

- the maximum steam sterilization temperature is 138°C / 280°F
- USA: the sterilization time (exposure time at sterilization temperature) is at least 4 minutes at least 132°C / 270°F.
- EU: the sterilization time (exposure time at sterilization temperature) is at least 4 minutes at least 134°C / 273°F.
- A minimum drying time of 30 minutes is recommended for each above cycles .

Care should be taken that:

- The steam sterilizers used meet the requirements of EN 13060 and EN 285.
- In general, use only adequately validated methods specific to the equipment and product for sterilization according to ISO 17865.
- Regularly check and service the sterilizers used.
- Observe the validated parameters in each cycle.

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