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SDS Swiss Dental Solutions AG  
Konstanzerstrasse 11  
CH-8260 Kreuzlingen

MOJE KERAMIK-IMPLANTATE GmbH & Co. KG  
Eschenweg 12  
07616 Petersberg

## Instructions for Use for Two-Piece SDS2.0 Dental Implants

**Description:** Two-piece SDS2.0 implants are inserted into the jawbone for transgingival healing. After implantation, the inner lumen of the two-piece SDS2.0 implant is covered with the enclosed cover screw made from the plastic PEEK (polyetheretherketone). For prosthodontic restoration after osseointegration the cover screw is removed, the abutment inserted, fixed with a screw and additionally cemented with glass ionomer cement (Ketac Cem™).

**Material:** SDS2.0 implants are made of high strength zirconia TZP-A (tetragonal zirconia polycrystal) zirconium dioxide ceramics (bio ceramics under ISO 13356).

**Indications SDS2.0:** Please check detailed clinical instructions in the SDS2.0 product information brochure. SDS2.0 dental implants are surgically inserted into the upper and lower jaw to replace missing teeth. They serve as attachment of prosthodontic appliances and are particularly suitable for patients with metal allergies and the chronic diseases resulting from them.

- Indicated for situations where implants are connected by interlocking or bridge restoration
- SDS1.1 4.6 / 5.4 mm Ø implants are permitted as single-tooth implant for incisivi, cuspids, premolars and molars and as bridge abutment
- SDS2.0 implants must be inserted at tissue-level

**Contraindications SDS2.0 implants:** Existing medical conditions or bad general state of health can limit the possibility to place dental implants surgically. Bruxism and insufficient bone quality/-quantity require specific measures to ensure treatment success. SDS implants are not suitable for applications in which the risk of excessive bending moments exists (e.g. extended crowns, extension bridges, bridges with more than one pontic unit).

1. SDS2.0 implants may not be used in bone-level position
2. SDS2.0 3.8 mm Ø implants may not be used as single tooth implants
3. SDS2.0 3.8 mm Ø implant may not be used in bridge restorations
4. Implant diameter smaller than 4.6 mm for upper central incisors, cuspids, molar region and/or bridge restoration
5. Bone not completely healed (residual otitis/NICO)
6. Serious and systemic health problems in patient
7. Bruxism
8. Untreated periodontitis, poor oral hygiene, untreated abscess or bone infection
9. Crown-length greater than the osseointegrated threaded segment
10. Cantilever bridges (mesial or distal)

11. Pontic width between two bridge abutments bigger than one pre-molar width
12. Connection of tooth with implant

**Warnings:** Anatomical and general health conditions can have a negative impact on dental implants. If the implant is exposed to pressure beyond its capacity, excessive bone loss or fracture of the implant can occur.

- The physician/dentist is obliged to inform patients about the following risk factors:
  - Poor bone quality
  - Poor oral hygiene
  - Diseases like blood disorders or untreated hormonal disorders
  - Alcohol or drug abuse
  - Two-piece SDS2.0 implants must be protected from load during healing

Observe the operating instructions of the unit manufacturer in laser applications.

**Possible adverse reactions:** The following adverse reactions are possible after insertion of a dental implant:


1. Pain, swelling, infection of soft- and hard tissue
2. Dysaesthesia/paraesthesia
3. No osseointegration
4. Loosening after osseointegration
5. Bone defects necessitating bone grafting
6. Perforation of the sinus, mandibular base, floor of the mouth or lower alveolar ducts
7. Damage of neighboring teeth/tooth roots
8. Excessive bone loss, which might necessitate surgical intervention
9. Esthetical problems
10. Fracture of the implant

**Precautionary measures:** Intensive diagnostics of the oral cavity must be performed ahead of every implantation. It is necessary to take and assess appropriate x-rays (OPG/DVT/CT) to clarify anatomical structures. When choosing the appropriate implant size/type, expected physiological chewing forces, as well as any parafunctional habits, must be taken into consideration. Special training is recommended to guarantee safe and effective insertion of ceramic SDS implants. Surgical standard procedures must be complied with, maximum insertion torque of 35 Ncm may not be exceeded. Incorrect surgical techniques can lead to functional failure of the implant and bone loss of the supporting bone structure or other side effects.

SDS2.0 implants may only be restored with SDS standard abutments. The maximum torque for the gold screw is 15 Ncm (PEEK: 5 Ncm). SDS implants may only be restored further with SDS products (or products which are authorized by SDS).

**Storage and handling:** SDS implants are supplied in sterile packaging and must be stored in their original packaging, under clean and dry conditions. If moisture is present in the inner package, or if the package is damaged or has been unintentionally opened, the implants must not be used. SDS implants must not be resterilized. Before use, the size/type and expiration date of the implants have to be checked against the label on the package. Appropriate hygienic measures must be observed when handling the implants. Contact with objects which could damage the implants must be avoided. SDS will not accept return of goods in damaged packages.

### GRAPHIC SYMBOLS:

	CATALOG NUMBER
	BATCH CODE
	CONSULT INSTRUCTIONS FOR USE
	STERILIZED USING STEAM
	DO NOT REUSE
	USE-BY DATE
	SYMBOL OF EUROPEAN CERTIFICATION NUM. 0483
	MANUFACTURER
	DO NOT RESTERILIZE
	TEMPERATURE LIMITATION

The detailed directions for use in the product brochure must be complied with without fail and are a basic requirement for the treatment. For technical support and further information please contact:

SDS Swiss Dental Solutions AG  
Konstanzerstrasse 11  
CH-8280 Kreuzlingen

MOJE KERAMIK-IMPLANTATE GmbH & Co. KG  
Eschenweg 12  
07616 Petersberg